

HISTONES (from Yeast cells)

The protocol below is given for the histone isolation from undegraded crude nuclear preparations and from whole cells, as developed from standard procedures in 1999 [see Waterborg (2000) J. Biol. Chem 275(17), 13007]. At 12 Mbp in the yeast genome, it is *assumed* that 0.012 pg/haploid genome yields, between G1 and G2 cells in log cultures, approx. 0.02 pg DNA (and thus histones) per nucleus or cell.

Histone purity is less than expected after isolation of nuclei, as judged by HPLC, SDS and AUT analysis. The level of H4 acetylation observed (Y1 preparation) is at least as high or up to 20% higher than the highest values reported. The best prep method reported is in Davie, Saunders, Walsh and Weber (1981) Nucleic Acids Res. 9: 3205-3216. It requires use of 100 mM butyrate during spheroplasting and nuclear isolation in all buffers for standard wild-type cells or at least 30 mM butyrate in low-protease (5%) strains. H4 Ac 0 to 4 = 14, 27, 31, 17 and 11 % (1.84 AcK/H4).

The observation of methionine oxidation of alfalfa H2B (see Binder 63) has prompted the inclusion of 2-mercaptoethanol in all buffers from the time that the cells are disrupted in the 40% GuCl homogenization step.

While it is possible to isolate intact histones from whole cells (see Y1-Y5 preps, 1999, binders 144 and 147), non-histone proteins that elute across a wide range of acetonitrile concentrations (30 to >50%) cause problems of contamination between samples, regenerating the column and bands in AU gels that obscure histone bands. Preparation Y5 established that, while HCl back-extraction of a whole cell extract is feasible and produces a fair yield (see details below), homogenization of cells by beading in a buffer selected to maximize nuclear and chromatin insolubility rather than nuclear purity, and minimize changes due to deacetylation or proteolysis, is effective. Acetylation levels were unchanged; no trace of proteolysis of histones. This method is described below as the Preferred method.

REQUIRED:

Nuclear Isolation Buffer (NIB)

The composition of this buffer has not been 'optimized' experimentally. It was based on the premise to create an isotonic, sucrose-based buffer that would stabilize nuclei and prevent solubilization of chromatin.

NIB composition: 0.25 M sucrose, 10 mM MgCl₂, 2.5 mM spermidine, 0.5 mM spermine, 20 mM HEPES, 100 mM butyrate, 0.1 % (w/v) Triton X-100, 5 mM 2-mercaptoethanol, pH 7.

NIB can be prepared 1 day before use and stored at 4°C until use.

PMSF to 1 mM is added just prior to use.

Required approx. 50(-100) ml per sample of up to 2 x 10¹⁰ cells; 500 ml for 8 samples (10¹⁰ cells)

Prepare: _____ ml NIB [100 ml] for _____ samples

_____ g sucrose [8.56 g]

_____ ml 1 M MgCl₂ stock [1 ml]

_____ µl spermidine stock [363 µl] (100 mg/ml in water, at MW=145.2 689 mM)

_____ µl spermine stock [101 µl] (100 mg/ml in water, at MW=202.3 494 mM)

_____ ml 1 M HEPES stock [2 ml]

_____ ml butyric acid [0.92 ml] (10.87 M solution)

_____ ml 25% Triton X-100 stock [0.40 ml]

_____ µl 2-mercaptoethanol (14.3 M) [35 µl]

Adjust pH to 7 by adding approx. 2.0-2.2 ml 5 N KOH per 100 ml.

Adjust with water to final volume. Storage (4°C): _____

Use of butyrate is based on Davie experience (NAR 9,3205,1981) despite the observation that butyrate added to live yeast cells does not affect acetylation level.

PMSF stock: freshly prepare 50 mM PMSF in isopropanol:

recommendation and cycloheximide, suggests that 10^{10} cells i.e. 200 ml of 5×10^7 cells/ml is required for H4 AUT gel patterns, with H3 then of all or 2/3 of H3 pool and with _ of H2B pool per AUT gel lane (long, large gels).

Cells in YPD (typically **YPD-S**) medium are used directly or concentrated from mid-log cultures (e.g. 5×10^7 cells/ml) by centrifugation e.g. 5 min $800 \times g$ (2000 rpm) at Room temperature. Use 250 ml conical Corning tubes, spin, decant and retain (or re-add) (used) YPD(-S) medium to (maximally) 10×10^7 cells/ml.

Notes on details of samples and/or incubations are below, or on separate pages.

Y6 (5/99): 1.25×10^{10} cells, incubated in Erlenmeijer flask, 30°C 150 rpm in incubator, with label in 90-100 ml YPD-S (pH approx 4) medium, grown for up to 2 hrs, thus likely up to 2.5×10^{10} cells, were collected as 2 plts in 50 ml tubes: plt volumes 1ml each, combined in NIB for beading.

2. Following incubation (or growth), collect cells from 50-(100 ml) samples into 1(or 2) marked 50 ml low-speed polypropylene tubes by centrifugation for 5 min at $800 \times g$ (2000 rpm). Alternatively, collect 200-250 ml into 1 Corning tube (250 ml), spin, resuspend in water and transfer to 50 ml tube for spin, decant, freeze and store at -80°C .
Decant and discard supernatant (or collect in radioactive waste container).

| | | | | | | | | | | | | |
|--------|--|--|--|--|--|--|--|--|--|--|--|--|
| ID | | | | | | | | | | | | |
| ml plt | | | | | | | | | | | | |

3. **Immediately** freeze the cell pellet (approx. 1 ml pellet from 10^{10} cells) in methanol-dry ice, and store (for at least 1 hour) at -80°C until isolation of histones.

Date: ____, ____-____-____, __:___ ____, ____-____-____, __:___

Isolation of (crude) nuclei (the “better” procedure to isolate histones)

4. Place 50 ml tubes with frozen cells on ice.
5. Use 3 ml **cold** Nuclear Isolation Buffer (NIB) to combine the cell pellets into one 50 ml tube with total sample volume of 4(-5) ml.
6. Add 500 micron glass beads (e.g. Sigma G9268, 425-600 μm , 30-40 U.S. sieve) until **just** no free liquid remains. This means that to 4 ml, glass is added to a total volume of approximately 10 ml. Homogenization on a fierce vortex of 10 ml in a 50 ml Corning tube is effective; a final volume of (15-)20 ml is TOO large and requires distribution across 2 tubes.
After the homogenate is collected (see below), the glass beads are extensively washed with water, rinsed on a Buchner filter, dried in an oven and collected for re-use.
7. Vortex for 1(-2) min vigorously the initially almost dry ball of glass beads which will start to liquify with vigorous semi-liquid eruptions of the surface.
8. Cool on ice.
9. To eliminate free standing liquid, add glass beads until no free liquid remains. The final volume may increase to approx. 13 ml.
10. Vortex for 1(-2) min vigorously. If initial vortex homogenization was effective, little or no increase in liquid state is observed
11. Add NIB to the 30 ml mark, mix and let the glass beads settle.
12. Decant the supernatant through 2 layers of Miracloth (in a funnel) into a Sorvall high-speed centrifugation tube (35 ml effective sample volume) **in ice**.
13. Add 10 ml NIB to the glass beads in the tube, mix and pour all into the funnel.
The glass beads are extensively washed with water, rinsed on a Buchner filter, dried in an oven and collected for re-use.

25. Add slowly with mixing: _____ μ l concentrated HCl (37%=12.1N) to a final concentration of 0.25 N.

The solution will likely turn slightly turbid. Observations:

26. Place the tubes in melting ice for (at least) 15 min to allow full precipitation of DNA and acidic, non-histone proteins: ____:____ - ____:____ hours.

Resulting color: _____

27. Clarify the solution by centrifugation for 30 min at 30,000(-40,000) x g (Sorvall SS34: 16,000 (to 18,000) rpm), 4°C and NO BRAKE.

cms note: Remember to put brake back on ☺

28. Carefully collect by decanting the clear supernatant from the small, white pellet.

(color: _____ size: _____)

If part of the precipitate, as pellet or floating on top of the solution, should end up in the supernatant, spin again to remove it before continuing to the next step.

29. Note the volume of the supernatant:

| | | | | | | | | | | | | | |
|-----|--|--|--|--|--|--|--|--|--|--|--|--|--|
| ID | | | | | | | | | | | | | |
| vol | | | | | | | | | | | | | |

30. Add 4 (-6) volumes of 0.1M KPi to decrease the GuCl concentration from (nominally) 40% to 5%.

| | | | | | | | | | | | | | |
|--------|--|--|--|--|--|--|--|--|--|--|--|--|--|
| ml KPi | | | | | | | | | | | | | |
|--------|--|--|--|--|--|--|--|--|--|--|--|--|--|

31. Measure the Refractive Index :

| | | | | | | | | | | | | | |
|----|--|--|--|--|--|--|--|--|--|--|--|--|--|
| n= | | | | | | | | | | | | | |
|----|--|--|--|--|--|--|--|--|--|--|--|--|--|

Add more 0.1M KPi (or 40% GuCl) to make the RI identical to that of 5% GuCl (~7.4 Brix).

| | | | | | | | | | | | | | |
|-----|--|--|--|--|--|--|--|--|--|--|--|--|--|
| KPi | | | | | | | | | | | | | |
| n= | | | | | | | | | | | | | |
| KPi | | | | | | | | | | | | | |
| n= | | | | | | | | | | | | | |
| KPi | | | | | | | | | | | | | |
| n= | | | | | | | | | | | | | |

32. Adjust the pH of the solution with 5N KOH to pH 6.8.

| | | | | | | | | | | | | | |
|-----|--|--|--|--|--|--|--|--|--|--|--|--|--|
| pH | | | | | | | | | | | | | |
| KOH | | | | | | | | | | | | | |

33. Add mercaptoethanol (!!) to 5 mM to this solution by adding 0.4 μ l/ml.

| | | | | | | | | | | | | | |
|--------------|--|--|--|--|--|--|--|--|--|--|--|--|--|
| Vol | | | | | | | | | | | | | |
| μ l 2-me | | | | | | | | | | | | | |

34. Skip this step if histones are extracted from nuclei.

In case of whole cell extraction, the solution has now become cloudy, especially using early to late stationary yeast cultures.

Clarify by centrifugation for 20 min at 5000 rpm at 4°C in Sorvall RC3 HG-4L or IEC PR6000 in conical 50 or 250 ml or large 1000 ml tubes.

Pellet observed: _____

43. The short 25 Gauge needle gives a correct, slow flow rate for column elution in a small volume.
44. Place a **collecting tube** below the syringe.
45. Elute the column with 10 column volumes (... ml) 40% GuCl.

| | | | | | | | | | | | | |
|----------------------------------|--|--|--|--|--|--|--|--|--|--|--|--|
| ID | | | | | | | | | | | | |
| # aliquots | | | | | | | | | | | | |
| ml 40% GuCl <i>vol eluate</i> | | | | | | | | | | | | |

46. Resin columns used with labeled nuclear preparations are discarded.
Columns with resin-only (if >5 ml resin volume) can be re-cycled/re-equilibrated.
47. Prepare Spectrapor 3 (MW cut-off 3500) for use by soaking in water.
18mm Spectrapor 3 holds 1 ml/cm; 45 mm Spectrapor 3 holds 6.4 ml/cm. Boiling in 1 mM EDTA for a few minutes is optional but not required and typically not done.
If low MW-cutoff membrane is unavailable, use 'normal' (~10000 MW cutoff) boiled in 5 mM EDTA in water for at least 10 minutes.
Cut pieces for 1 or 2 samples, bordered by knots (5 cm/knot), 1 knot between 2 samples.
48. Pipet or pour (using a funnel) the eluate into the dialysis bag.
- Close each end by a double knot.
 - Do NOT exclude all air. It will assure a vertical position of all dialysis bags during dialysis, perpendicular to the stirred current.
 - Use simple rubber bands to tie each bag to a marked, identified (inverted) pipet.
Multiple sets of small samples, up to 2 per tubing, separated by a single knot, with clearly marked order on the pipet, may be tied to a single holder.
49. Initial dialysis: 1 hour against 100 volumes of 2.5% (v/v) acetic acid (with approx. 1 mM 2-mercaptoethanol: 0.1 µl/ml or 100 µl/liter) at 4°C, to remove the major amount of GuCl.
The use of 0.5, 1.0, 2.0 or 5.0 liter polypropylene beakers with handle, with large magnetic flea, on magnetic stirrer, in a 4°C cold-room works very well.
_____ ml 2.5% HAc/2me. Date: ____, ___-___-____:___ h. ____, ___-___-____:___ h.
50. Dialysis 2: 1 hour against 100 volumes.
_____ ml 2.5% HAc/2me. Date: ____, ___-___-____:___ h. ____, ___-___-____:___ h.
51. Dialysis 3 : overnight against 100 volumes.
_____ ml 2.5% HAc/2me. Date: ____, ___-___-____:___ h. ____, ___-___-____:___ h.
- Ultra-filtration through Centricon-100 (2 ml sample volume) or CentriPlus-100 (15 ml sample volume) through Amicon YM-100 membrane (100,000 D MW cut-off) gives a filtrate with quantitative histone yield, removing “insolubles” and compounds that elute heterogeneously at high acetonitrile concentrations (30% - >60%) in reversed-phase ProteinPlus HPLC.
Incorporation of total label into dialysate is reduced to 25% (*Ac 2h), 3% (*Ac 2h +cycloheximide), 50-60% (*arg, lys 2h), 15-20% (*met 2hr ±cycloheximide), 10% (*AdoMet 2hr ±cycloheximide) (Y9 experiment), apparently retaining all label incorporated in yeast histones.

52. Transfer the dialysate(s) to CODED/labeled, capped polypropylene tube(s) (15, 50 or 250 ml).

| | | | | | | | | | | | | |
|------------|--|--|--|--|--|--|--|--|--|--|--|--|
| ID | | | | | | | | | | | | |
| total (ml) | | | | | | | | | | | | |
| for LSC: | | | | | | | | | | | | |

53. Centricon-100 filtration of dialysates in aliquots of up to 2 ml at 4°C and 1000xg (2200 rpm) in rubber adapter-fitted IEC PR6000 (or RC-2C HG-4L rotor with adapters).

- initial 2 ml spin 60 min: _____ - _____ h. Result: _____
- load more/rest of dialysate: _____ - _____ h. Result: _____
- further spin(s): _____ - _____ h. Result: _____
- ONE wash of the small residual volume with (up to) 2 ml 2.5% HAc and spin to “completion” (overnight). Date: ____, ___-___-___ __:___ h. ____, ___-___-___ __:___ h.

54. Transfer the filtrate(s) to CODED/labeled, capped polypropylene tube(s) (15, 50 or 250 ml).

| | | | | | | | | | | | | |
|------------|--|--|--|--|--|--|--|--|--|--|--|--|
| ID | | | | | | | | | | | | |
| total (ml) | | | | | | | | | | | | |
| for | | | | | | | | | | | | |
| for | | | | | | | | | | | | |
| for | | | | | | | | | | | | |
| for | | | | | | | | | | | | |

Aliquots of these dialysates may be taken for lyophilization in separate tubes, e.g. in 1.5 or 0.5 ml Eppendorf tubes, for uses as sample in single HPLC or gel analyses.

55. Puncture the caps of the tubes.

56. Freeze the tubes in -70°C freezer (or with methanol-dry ice or liquid nitrogen.)

Date: ____, ___-___-___ __:___ hours. ____, ___-___-___ __:___ hours.

Do NOT place tube in insulating (polystyrene) holders: this causes melting and prevent sublimation.

57. Lyophilize until dry in an oil-pump vacuum with methanol-filled cold-trap and chemical trap.

Despite all traps, do NOT use the oil pump for more than 250 hours (after heavy use, especially drying gels, less !) without **completely** changing the oil (disassemble and wipe all accessible parts in the Savant oil pump). (Do NOT use the so-called ‘cleaning solvent’: it causes the pump to seize!)

Lyophilization. Date: ____, ___-___-___ __:___ hours. ____, ___-___-___ __:___ hours.

Result: _____ Continue ____, ___-___-___ __:___ hours.

Result: _____ Continue ____, ___-___-___ __:___ hours.

Observations:

Dried histone lyophilizates can be stored indefinitely at -20°C in closed (new caps!) tubes, especially under nitrogen atmosphere. Even storage at room temperature is acceptable.

NOTE : the residual mercaptoethanol in the lyophilizate will confer protection against oxidation if excessive air contact is avoided: store dry histones in tightly capped tubes !